

User-Developed Protocol:

Purification of total RNA from fatty tissues using QIAzol Lysis Reagent, MaXtract High Density, and the TissueRuptor

This protocol has been adapted by customers for the purification of total RNA from fatty tissues using QIAzol Lysis Reagent, MaXtract High Density, and the TissueRuptor. This protocol can also be used with all other types of tissues. **This protocol has not been thoroughly tested and optimized by QIAGEN.**

IMPORTANT: Please be sure to read the *QIAzol Handbook* and the *MaXtract Low and High Density Handbook* before starting, paying careful attention to the Safety Information section.

Equipment and reagents to be supplied by user

- QIAzol Lysis Reagent (200 ml) (cat. no. 79306)
- MaXtract High Density (200 x 2 ml) (cat. no. 129056)
- TissueRuptor (to order, visit www.qiagen.com or contact your local QIAGEN office)
- Also refer to “Equipment and Reagents to Be Supplied by User” in the *QIAzol Handbook*

Important points before starting

- Ensure that you are familiar with operating the TissueRuptor by referring to the *TissueRuptor User Manual* and *TissueRuptor Handbook*. For other rotor–stator homogenizers, refer to suppliers’ guidelines.
- If using QIAzol Lysis Reagent for the first time, read “Important Notes” in the *QIAzol Handbook* (October 2006 edition).
- Fresh, frozen, or RNA^{later}[®] stabilized tissues can be used.* If freezing tissues, flash-freeze in liquid nitrogen and immediately transfer to –70°C, where they can be stored for several months. Do not allow tissues to thaw during weighing or handling prior to disruption in QIAzol Lysis Reagent. Homogenized tissue lysates from step 3 can also be stored at –70°C for at least 1 month. Incubate frozen lysates at 37°C in a water bath until completely thawed and salts are dissolved before continuing with step 4. Avoid prolonged incubation, which may compromise RNA integrity.

Procedure

- 1. Add QIAzol Lysis Reagent to an appropriate vessel for disruption and homogenization and subsequent centrifugation: 1 ml QIAzol Lysis Reagent per 100 mg tissue is required. The volume of tissue should not exceed 10% of the volume of QIAzol Lysis Reagent.**

Generally, round-bottomed tubes allow more efficient disruption and homogenization than conical-bottomed tubes.

* RNA^{later} RNA Stabilization Reagent cannot be used with adipose tissue due to the high abundance of fat, but can be used with other fatty tissues such as brain.

- 2. Excise the tissue sample from the animal or remove it from storage. Determine the amount of tissue and place it into the QIAzol Lysis Reagent. Proceed immediately to step 3.**

Weighing tissue is the most accurate way to determine the amount.

If the tissue sample was stored in RNA*later* RNA Stabilization Reagent, remove it from the reagent using forceps and be sure to remove any crystals that may have formed.

RNA in harvested tissues is not protected until the tissues are treated with RNA*later* RNA Stabilization Reagent, flash-frozen, or disrupted and homogenized in step 3. Frozen tissues should not be allowed to thaw during handling. The relevant procedures should be carried out as quickly as possible.

- 3. Place the tip of the TissueRuptor disposable probe into the QIAzol Lysis Reagent, and operate the TissueRuptor at full speed until the tissue lysate is uniformly homogeneous (usually 20–40 s).**

Note: To avoid damage to the TissueRuptor and disposable probe during operation, make sure the tip of the probe remains submerged in the QIAzol Lysis Reagent.

Note: Incomplete homogenization leads to significantly reduced RNA yields. Homogenization with the TissueRuptor or TissueLyser generally results in higher RNA yields than with other methods.

Foaming may occur during homogenization, especially of brain tissue. If this occurs, let the homogenate stand at room temperature for 2–3 min until the foam subsides before continuing with the procedure.

Optional: For samples containing a relatively high content of fat, proteins, polysaccharides, or extracellular material, centrifuge the homogenate at 12,000 x g for 10 min at 4°C to remove insoluble material. Carefully transfer the supernatant to a new tube, and proceed to step 4.

- 4. Place the tube containing the homogenate on the benchtop at room temperature (15–25°C) for 5 min.**

This step promotes dissociation of nucleoprotein complexes.

- 5. Centrifuge a 2 ml tube of MaXtract High Density at 12,000–16,000 x g for 20–30 s. Transfer the homogenate to the MaXtract tube.**

Note: If the homogenate was derived from RNA*later* stabilized tissue, add 50–100 µl RNase-free water to the MaXtract tube as well.

- 6. To the MaXtract tube, add 0.2 ml chloroform per 1 ml QIAzol Lysis Reagent pipetted in step 1. Securely cap the tube, and shake it vigorously for 15 s. Do not vortex.**

Thorough mixing is important for subsequent phase separation.

- 7. Place the MaXtract tube on the benchtop at room temperature for 2–3 min.**

- 8. Centrifuge at 12,000 x g for 15 min at 4°C.**

After centrifugation, the sample separates into 3 phases: an upper, colorless, aqueous phase containing RNA; a white interphase (MaXtract High Density gel); and a lower, red, organic phase. For tissues with an especially high fat content, an additional clear phase may be visible below the red, organic phase. The volume of the aqueous phase is approximately 60% of the volume of the QIAzol Lysis Reagent pipetted in step 1.

- 9. Transfer the upper, aqueous phase to a new tube. Add 0.5 ml isopropanol per 1 ml QIAzol Lysis Reagent pipetted in step 1. Mix thoroughly by vortexing.**

- 10. Place the tube on the benchtop at room temperature for 10 min.**
- 11. Centrifuge at 12,000 x g for 10 min at 4°C.**
- 12. Carefully aspirate and discard the supernatant.**

The RNA pellet is often visible as a gel-like or white pellet at the bottom of the tube.
- 13. Add at least 1 ml of 75% ethanol per 1 ml QIAzol Lysis Reagent pipetted in step 1. Centrifuge at 7500 x g for 5 min at 4°C.**

If the RNA pellet floats or sticks to the side of the tube, bring it to the bottom of the tube by centrifuging at 12,000 x g for 5 min at 4°C.
- 14. Remove the supernatant completely, and briefly air-dry the RNA pellet.**

Do not dry the RNA using a vacuum.
- 15. Redissolve the RNA in an appropriate volume of RNase-free water. Clean up the RNA using the RNeasy[®] MinElute[®] Cleanup Kit or RNeasy Mini, Midi, or Maxi Kit.**

We recommend RNA cleanup to remove contaminating phenol. The RNeasy MinElute Cleanup Kit and RNeasy Mini, Midi, and Maxi Kits allow cleanup of up to 45 µg, 100 µg, 1 mg, and 6 mg total RNA, respectively. For details, refer to the RNA cleanup protocol in the handbook supplied with these kits.

QIAzol Lysis Reagent, RNA^{later} RNA Stabilization Reagent, the RNeasy MinElute Cleanup Kit, the RNeasy Midi Kit, and the RNeasy Maxi Kit are intended for research use. No claim or representation is intended to provide information for the diagnosis, prevention, or treatment of a disease. The RNeasy Mini Kit is intended for general laboratory use. No claim or representation is intended to provide information for the diagnosis, prevention, or treatment of a disease.

QIAGEN[®] kit handbooks can be requested from QIAGEN Technical Service or your local QIAGEN distributor. Selected kit handbooks can be downloaded from www.qiagen.com/literature. Material safety data sheets (MSDS) for any QIAGEN product can be downloaded from www.qiagen.com/ts/msds.asp.

Trademarks: QIAGEN[®], MinElute[®], RNeasy[®] (QIAGEN Group).

“RNA^{later}” is a trademark of AMBION, Inc., Austin, Texas and is covered by various U.S. and foreign patents.

QIAzol Lysis Reagent is a subject of US Patent No. 5,346,994 and foreign equivalents.

© 2007 QIAGEN, all rights reserved.

Dear Customer,

Thank you for downloading this supplementary protocol. We hope it is helpful for you.

In order to improve the quality of our service, we would appreciate your feedback.

Please submit your comments or suggestions at:

www.qiagen.com/ProtocolFeedback

Thank you for choosing QIAGEN

Best regards,

Your Technical Service Team

