

NanoAmp™ RT-IVT Labeling Kit

RT-IVT Labeling Process for **One Round** of Amplification

Safety

For safety and biohazard guidelines, refer to the “Safety” section in the *NanoAmp™ RT-IVT Labeling Kit Protocol* (PN 4365710). For all chemicals that appear in **bold** type throughout this document, read the MSDS and follow the handling instructions. Wear appropriate protective eyewear, clothing, and gloves.


Overview

The purpose of the NanoAmp™ RT-IVT Labeling Kit (PN 4365715) is to convert RNA into digoxigenin (DIG)-labeled cRNA for hybridization to Applied Biosystems microarrays. You can use the kit to perform:

- One round of amplification (standard)
or
- Two rounds of amplification

This Quick Reference Card briefly describes how to use the kit for one round of amplification. For two rounds of amplification, refer to Quick Reference Card PN 4370891. For a detailed protocol of both processes, refer to the *NanoAmp™ RT-IVT Labeling Kit Protocol* (PN 4365710).

Workflow for One Round of Amplification



Estimated completion time: 5–6 hours plus one overnight incubation.

Handling the Enzymes

IMPORTANT! Keep all enzymes on ice at all times throughout these procedures. Do not vortex the enzymes.

Performing Reverse Transcription

Required Materials

- RNA sample: up to 2 µg of total RNA (for one round of amplification)
- Reverse transcription reagents (blue-capped tubes):
 - Control RNA
 - T7-Oligo (dT) Primer
 - **10× 1st Strand Buffer**
 - dNTP Mix
 - **RT Enzyme**
 - **RNase Inhibitor**
- Nuclease-free water (gray-capped tube)
- 0.2-mL MicroAmp® reaction tube(s)
- Ice bucket
- Pipettor: 1- to 20-µL range
- Pipette tips
- Vortexer
- Microcentrifuge
- Thermal cycler

Procedure

1. Dilute the control RNA 1:5000:
 - a. Add 2 µL of control RNA to 98 µL of nuclease-free water, then mix by vortexing (1:50).
 - b. Add 2 µL of the first dilution to 198 µL of nuclease-free water, then mix by vortexing (1:100, final is 1:5000).
2. Thaw the RNA sample and all reagents on ice.
3. Vortex separately the RNA sample, T7-Oligo (dT) Primer, **10× 1st Strand Buffer**, and dNTP Mix.
4. Centrifuge all reagent tubes briefly.
5. In the order given below, add the components to a 0.2-mL MicroAmp reaction tube on ice, then mix by pipetting up and down:

Component	Volume (µL)
T7-Oligo (dT) Primer	1.0
Diluted control RNA (1:5000)	2.0
RNA sample (up to 2 µg of total RNA); bring volume to 9.0 µL with nuclease-free water	9.0
Total	12.0

6. Heat and then cool the mixture in a thermal cycler:

Stage	Temp. (°C)	Time	Reaction Vol.
1	70	5 min	12 µL
2	4	Indefinite hold	

7. After the run, place the reaction tube on ice.
8. If you see precipitate in the **10× 1st Strand Buffer**, warm it at 37 °C for 5 min, then vortex briefly.
9. In the order given below, add the components to the reaction tube on ice, then mix thoroughly by pipetting up and down:

Component	Volume (µL)
T7-Oligo (dT) Primer, diluted control RNA, and RNA sample	12.0
10× 1st Strand Buffer	2.0
dNTP Mix	4.0
RT Enzyme	1.0
RNase Inhibitor	1.0
Total	20.0

Note: Alternatively, you can prepare a master mix of the **10× 1st Strand Buffer**, dNTP Mix, **RT Enzyme**, and **RNase Inhibitor** up to 15 min before use. Add the primer/RNA mixture to an 8-µL aliquot of this master mix. Add sufficient overage for all reactions. Mix gently, then go to [step 10](#).

10. Perform reverse transcription in a thermal cycler:

Stage	Temp. (°C)	Time	Reaction Vol.
1	25	10 min	20 µL
2	42	2 hours	
3	70	5 min	
4	4	Indefinite hold	

11. After the run, place the reaction tube on ice.

Performing Second-Strand Synthesis

Required Materials

- Second-strand synthesis reagents (yellow-capped tubes):
 - **10× 2nd Strand Buffer**
 - **DNA Polymerase**
 - **RNase H**
- dNTP Mix (blue-capped tube)
- Nuclease-free water (gray-capped tube)
- Ice bucket
- Pipettors: 1- to 20-µL range, 20- to 200-µL range
- Pipette tips
- Thermal cycler

Procedure

- Thaw all reagents on ice.
- Vortex separately the **10× 2nd Strand Buffer** and dNTP Mix.

3. Centrifuge all reagent tubes briefly.

- If you see precipitate in the **10× 2nd Strand Buffer**, warm it at 37 °C for 5 min, then vortex briefly.
- In the order given below, add the components to the cDNA mixture on ice, then mix gently:

Component	Volume (µL)
cDNA mixture	20.0
Nuclease-free water	63.0
10× 2nd Strand Buffer	10.0
dNTP Mix	4.0
DNA Polymerase	2.0
RNase H	1.0
Total	100.0

Note: Alternatively, you can prepare a master mix of the water, **10× 2nd Strand Buffer**, dNTP Mix, **DNA Polymerase**, and **RNase H** up to 15 min before use. Add the cDNA mixture to an 80-µL aliquot of this master mix. Add sufficient overage for all reactions. Mix gently, then go to [step 6](#).

6. Perform second-strand synthesis in a thermal cycler:

Stage	Temp. (°C)	Time	Reaction Vol.
1	16	2 hours	100 µL
2	70	5 min	
3	4	Indefinite hold	

7. After the run, place the reaction tube on ice.

Purifying cDNA

Required Materials

- cDNA purification components:
 - **DNA Binding Buffer**
 - **Wash Buffer**
 - Nuclease-free water
 - DNA purification column(s) and tube(s)
 - DNA elution tube(s)
- 100% ethanol**, ACS-grade or better
- Pipettors: 20- to 200-µL range, 100- to 1000-µL range
- Pipette tips
- Microcentrifuge

IMPORTANT! Before you use the kit for the first time, reconstitute the **Wash Buffer** by adding 16 mL of **100% ethanol** to the **Wash Buffer** bottle. Mix well, then mark the label to indicate that **100% ethanol** was added.

Procedure

- If you see precipitate in the **DNA Binding Buffer**, warm it at 37 °C for 5 min, then vortex briefly.
- In a new 1.5-mL nuclease-free microcentrifuge tube:
 - Combine the **DNA Binding Buffer** (250 µL) and the entire second-strand synthesis reaction (100 µL).

b. Mix thoroughly by pipetting up and down.

IMPORTANT! Do not use the DNA elution tubes provided in the kit until [step 6](#). Use other nuclease-free microcentrifuge tubes for mixing the **DNA Binding Buffer** and second-strand synthesis reaction.

3. Begin the cDNA purification:

IMPORTANT! Before using the columns, inspect each one to be sure the filter is flush against the bottom of the column. If necessary, push it into place with a pipette tip.

- Seat the DNA purification column in its wash tube.
- Pipette the **DNA Binding Buffer**/reaction mixture (350 μ L) onto the DNA purification column.
- Centrifuge the column and tube at 10,000 $\times g$ for 1 min.
- Make sure that the entire volume passes through the column. If it does not, centrifuge the column and tube at 10,000 $\times g$ for 1 min.
- Remove the column from the tube, discard the liquid, then reinsert the column into the tube.

4. Wash the cDNA:

- Add 500 μ L of **Wash Buffer** (containing **100% ethanol**) to the column.
- Centrifuge the column and tube at 10,000 $\times g$ for 1 min.
- Remove the column from the tube, discard the liquid, then reinsert the column into the tube.

5. Centrifuge the column and tube again at 10,000 $\times g$ for 1 min.

6. Elute the cDNA (elution #1):

- Transfer the column to a new DNA elution tube ([use the tubes provided in the kit](#)).
- Pipette 10 μ L of nuclease-free water onto the center of the fiber matrix at the bottom of the column.
- Incubate the column at room temperature for 2 min.
- Centrifuge the column and tube at 10,000 $\times g$ for 1 min for an elution volume of approximately 9 μ L.

7. Elute the cDNA (elution #2):

- Pipette 10 μ L of nuclease-free water onto the center of the fiber matrix at the bottom of the column.
- Incubate the column at room temperature for 2 min.
- Centrifuge the column and tube at 10,000 $\times g$ for 1 min for a total elution volume of approximately 18 μ L.

8. Discard the column, then close the tube.

Note: If desired, you can stop at this point and store the purified cDNA at $-20\text{ }^{\circ}\text{C}$ overnight.

Performing IVT Labeling

Required Materials

- IVT labeling reagents (green-capped tubes):
 - 10 \times IVT Buffer**
 - NTP Mix**
 - IVT Control DNA
 - IVT Enzyme Mix**
- DIG-UTP (You need two tubes of DIG-UTP for one labeling kit [Roche Molecular Biochemicals, Cat. No. 03359247910].)
- Nuclease-free water (gray-capped tube)

- Pipettors: 1- to 20- μ L range, 20- to 200- μ L range
- Pipette tips
- Thermal cycler

Procedure

- Thaw the **10 \times IVT Buffer** at room temperature; thaw the remaining reagents on ice.
- Vortex separately the **10 \times IVT Buffer**, **NTP Mix**, IVT Control DNA, and DIG-UTP.
- Centrifuge all reagent tubes briefly.
- If you see precipitate in the **10 \times IVT Buffer**, warm the buffer at 37 $^{\circ}\text{C}$ for 5 min, then vortex briefly.
- In the order given below, add the components to the cDNA output at room temperature. Mix gently, then briefly centrifuge the tube.

Component	Volume (μ L)
ds cDNA output; bring volume to 18 μ L with nuclease-free water	18.0
10\times IVT Buffer	4.0
DIG-UTP	8.0
NTP Mix	4.0
IVT Control DNA	2.0
IVT Enzyme Mix	4.0
Total	40.0

Note: Alternatively, you can prepare a master mix of the **10 \times IVT Buffer**, DIG-UTP, **NTP Mix**, IVT Control DNA, and **IVT Enzyme Mix** up to 15 min before use. Add the ds cDNA output to a 22- μ L aliquot of this master mix. Add sufficient overage for all reactions. Mix gently, then go to [step 6](#).

6. Perform IVT in a thermal cycler:

Stage	Temp. ($^{\circ}\text{C}$)	Time	Reaction Vol.
1	37	9 hours	40 μ L
2	4	Indefinite hold	

7. After the run, remove the tube from the thermal cycler.

Purifying cRNA

Required Materials

- cRNA purification components
 - RNA Binding Buffer**
 - Wash Buffer**
 - Nuclease-free water
 - RNA purification column(s)
 - RNA collection tube(s)
- 100% ethanol**, ACS-grade or better
- Ice bucket
- 1.5-mL nuclease-free microcentrifuge tube(s)
- Pipettors: 20- to 200- μ L range, 100- to 1000- μ L range
- Pipette tips

- Microcentrifuge

Procedure

1. In a new 1.5-mL nuclease-free microcentrifuge tube:
 - a. Combine the nuclease-free water (60 μ L) and the entire IVT reaction (40 μ L).
 - b. Vortex briefly to mix.

IMPORTANT! Do not use the RNA collection tubes provided in the kit until [step 6](#). Use other nuclease-free microcentrifuge tubes for mixing the water and IVT reaction.

2. To the water/IVT reaction mixture:
 - a. Add the **RNA Binding Buffer** (350 μ L) and the **100% ethanol** (250 μ L).
 - b. Mix by pipetting up and down (do not vortex).

3. Begin the cRNA purification:

IMPORTANT! Before using the columns, inspect each one to be sure the filter is flush against the bottom of the column. If necessary, push it into place with a pipette tip.

- a. Insert an RNA purification column into an RNA collection tube.
 - b. Pipette the IVT reaction/**RNA Binding Buffer**/**ethanol** mixture (700 μ L) onto the column, then close the tube.
 - c. Centrifuge the column and tube at $10,000 \times g$ for 1 min.
 - d. Make sure that the entire volume passes through the column. If it does not, centrifuge the column and tube at $10,000 \times g$ for 1 min.
 - e. Remove the column from the tube, discard the liquid, then reinsert the column into the tube.
4. Wash the cRNA:
 - a. Add 650 μ L of **Wash Buffer** to the column, then close the tube.
 - b. Centrifuge the column and tube at $10,000 \times g$ for 1 min.
 - c. Remove the column from the tube, discard the liquid, then reinsert the column into the tube.
 5. Close the tube, then centrifuge the column and tube at $10,000 \times g$ for 1 min.
 6. Elute the cRNA:
 - a. Transfer the column to a new RNA collection tube ([use the tubes provided in the kit](#)).
 - b. Pipette 100 μ L of nuclease-free water onto the fiber matrix at the bottom of the column, then close the tube.
 - c. Incubate the column at room temperature for 2 min.
 - d. Centrifuge the column and tube at $10,000 \times g$ for 1 min.
 7. Discard the column, then close the tube containing the eluted RNA.
 8. Store the cRNA product on ice while you assess its quantity and quality.

Assessing cRNA Quantity and Quality

Assess the quantity and quality of the cRNA product before proceeding with the Applied Biosystems Chemiluminescence Detection Kit for Applied Biosystems gene expression microarrays.

Required Materials

- TE buffer (10 mM Tris, 1 mM EDTA, pH 8.0), user-supplied
- Pipettors: 1- to 20- μ L range, 20- to 200- μ L range
- Pipette tips
- Cuvette
- UV spectrophotometer
- 1 to 2% agarose gel

Procedure

1. Measure the ultraviolet absorbance:
 - a. Gently mix the cRNA sample.
 - b. Dilute a small amount of cRNA product 1:30 with TE buffer.
 - c. Measure the absorbance at 260 nm.
2. Calculate the concentration and yield:
 - cRNA concentration (μ g/ μ L) = $A_{260} \times 0.04 \mu\text{g}/\mu\text{L} \times 30$ (dilution factor)
 - cRNA yield (μ g) = cRNA conc (μ g/ μ L) \times 100 μ L total volume of cRNA purified
3. Analyze the cRNA: Prepare a 1 to 2% agarose gel, then run 0.5 to 1 μ g of cRNA through the gel to visualize the cRNA product.

Note: You can also analyze an aliquot of cRNA with the Agilent 2100 Bioanalyzer. Refer to the manufacturer's instructions for further information.

(Optional) Concentrating the Purified cRNA

If necessary, concentrate the cRNA by vacuum centrifugation. Use only the medium- or low-temperature settings during the drying process. Remove the sample when the desired volume is achieved.

Storing cRNA Product

- -15 to -25 $^{\circ}\text{C}$ for up to 2 months
- -80 $^{\circ}\text{C}$ for long-term storage

© Copyright 2006, Applied Biosystems. All rights reserved. For Research Use Only. Not for use in diagnostic procedures.

NOTICE TO PURCHASER: Please refer to the *NanoAmp™ RT-IVT Labeling Kit Protocol* for limited label license or disclaimer information.

TRADEMARKS: AB (Design), Applied Biosystems, and MicroAmp are registered trademarks and Applera and NanoAmp are trademarks of Applera Corporation and/or its subsidiaries in the U.S. and/or certain other countries.

All other trademarks are the sole property of their respective owners.

Printed January 2006

